

TGFBR1 antibody



Catalog Number: 116023

Product name

TGFBR1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TGFBR1 Rabbit Polyclonal antibody. Positive WB detected in A549 cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004612).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116023(TGFBR1 Antibody) at dilution of 1:300