

TMEM98 antibody



Catalog Number: 116200

Product name

TMEM98 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TMEM98 Rabbit Polyclonal antibody. Positive WB detected in mouse E12.5 embryos, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 21-25 kDa

Preparation

This antibody was obtained by immunization of TMEM98 recombinant protein (Accession Number: NM_015544). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: N/A

Validations



WB result of anti-TMEM98(Catalog No:116200) in WT and KO mouse (10 ug of total protein lysate from E12.5 embryos) by Dr. Sally H. Cross.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116200(TMEM98 Antibody) at dilution of 1:1000