

UBE4A antibody



Catalog Number: 116543

Product name

UBE4A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

UBE4A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse brain tissue, L02 cells. Observed molecular weight by Western-blot: 123-130kd

Preparation

This antibody was obtained by immunization of UBE4A recombinant protein (Accession Number: NM_001204077). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116543(UBE4A antibody) at dilution of 1:500