

Anti-USP13 Rabbit antibody



Catalog Number: 168826

Product name

Anti-USP13 Rabbit antibody

Specificity

Human, Mouse

Antibody description

USP13 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human USP13

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

FC: 1/100-1/500

IP: 1/20-1/50

Validations



Western blot detection of USP13 in U2OS cell lysates using USP13 Rabbit mAb(1:500 diluted).Predicted band size:97kDa.Observed band size:97kDa.