

Anti-CD64/FCGR1 antibody



Catalog Number: 100179

Product name

Anti-CD64/FCGR1 antibody

Immunogen

[Rat CD64/FCGR1 \(His Tag\) recombinant protein](#)

Specificity

Rat CD64 / FCGR1

Antibody description

Rabbit monoclonal to CD64/FCGR1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat CD64 / FCGR1 (rR CD64 / FCGR1; NP_001094306.1; Met1-Pro285) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD64 / FCGR1 Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Rat FCGR1(CD64) expression on SD rat bone marrow cells. Cells were stained with FITC-conjugated anti-Rat FCGR1(CD64). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.