

Anti-RAB9A / RAB9 antibody



Catalog Number: 100268

Product name

Anti-RAB9A / RAB9 antibody

Specificity

Human RAB9A

Antibody description

Rabbit polyclonal to RAB9A / RAB9

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human RAB9A, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 10-20 µg/ml

Validations



Lanes	A	B	C
Items			
Sample (whole cell lysate)	Jurkat	K562	HEK293
Sample Volume (µg/lane)	30	30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	10-20 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Explanation	Developed using Odyssey imaging system. Predicted band size : 28 kDa Observed band size : 28 kDa Additional bands at : 65 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).		

RAB9A / RAB9 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot