

Anti-PTH / PTH1 / Parathyroid Hormone antibody



Catalog Number: 100280

Product name

Polyclonal

Anti-PTH / PTH1 / Parathyroid Hormone antibody

Ig Type

Rabbit IgG

Immunogen

[Human PTH / PTH1 / Parathyroid Hormone \(aa 32-65, GST Tag\) recombinant protein](#)

Applications

IHC-P

Specificity

Human PTH / PTH1 / Parathyroid Hormone

Dilutions

IHC-P: 1-10 µg/ml


Antibody description

Rabbit polyclonal to PTH / PTH1 / Parathyroid Hormone

Validations

Preparation

Produced in rabbits immunized with purified, recombinant Human PTH / PTH1 / Parathyroid Hormone (rh PTH / PTH1 / Parathyroid Hormone; P05062; Ser32-Phe65). PTH / PTH1 / Parathyroid Hormone specific IgG was purified by Human PTH / PTH1 / Parathyroid Hormone affinity chromatography.



Formulation

0.2 µm filtered solution in PBS


PTH / PTH1 / Parathyroid Hormone Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.



Clonality

PTH / PTH1 / Parathyroid Hormone Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry