

Anti-RACK1-GNB2L1 antibody



Catalog Number: 100324

Product name

Anti-RACK1-GNB2L1 antibody

Immunogen

[Human RACK1-GNB2L1 \(His & MBP Tag\) recombinant protein](#)

Specificity

Human RACK1 / GNB2L1

Antibody description

Rabbit polyclonal to RACK1-GNB2L1

Preparation

Produced in rabbits immunized with purified, recombinant Human RACK1 / GNB2L1 (rh RACK1 / GNB2L1; P63244; Met1-Arg317). RACK1 / GNB2L1 specific IgG was purified by Human RACK1 / GNB2L1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human RACK1 / GNB2L1. The detection limit for Human RACK1 / GNB2L1 is approximately 0.00975 ng/well.

IHC-P: 0.1-2 µg/mL

Validations



RACK1 / GNB2L1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human GNB2L1 in human hepatoma (from 2 donors) with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).