

# Anti-Epcr/PROCR antibody



Catalog Number: 100464

## Product name

Anti-Epcr/PROCR antibody

## Immunogen

[Cynomolgus Epcr/PROCR \(His Tag\) recombinant protein](#)

## Specificity

Cynomolgus Epcr / PROCR

## Antibody description

Rabbit polyclonal to Epcr/PROCR

## Preparation

Produced in rabbits immunized with purified, recombinant Cynomolgus Epcr / PROCR (Met1-Thr209). Epcr / PROCR specific IgG was purified by Cynomolgus Epcr / PROCR affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 0.1-2 µg/mL

## Validations



Epcr / PROCR Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of cynomolgus PROCR in cynomolgus kidney with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).



Epcr / PROCR Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of cynomolgus PROCR in cynomolgus lung with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).