

# Anti-BLyS/TNFSF13B antibody



Catalog Number: 100510

## Product name

Anti-BLyS/TNFSF13B antibody

## Immunogen

[Human BLyS/TNFSF13B recombinant protein](#)

## Specificity

Human BLyS / TNFSF13B / BAFF

## Antibody description

Rabbit monoclonal to BLyS/TNFSF13B

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human BLyS / TNFSF13B / BAFF (rh BLyS / TNFSF13B / BAFF; Q9Y275-1; Ala134-Leu285).

## Formulation

0.2 µm filtered solution in PBS with 5% treha

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions


WB: 5-10 µg/ml

IP: 1-4 µg/mg of lysate

## Validations



BLyS / TNFSF13B / BAFF Antibody, Rabbit MAb,



Western blot

BLyS / TNFSF13B / BAFF Antibody, Rabbit MAb,  
Immunoprecipitation