

# Anti-RBP4 antibody



Catalog Number: 100553

## Product name

Anti-RBP4 antibody

## Immunogen

[Mouse RBP4 \(His Tag\) recombinant protein](#)

## Specificity

Mouse RBP4

## Antibody description

Rabbit monoclonal to RBP4

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse RBP4 (rM RBP4; NP\_035385.1; Met1-Leu201).

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IF, ICC/IF

## Dilutions

ICC/IF: 10-25  $\mu$ g/mL

## Validations



RBP4 Antibody, Rabbit MAb, Immunofluorescence

Immunofluorescence staining of mouse RBP4 in mouse splenocytes. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with rabbit anti-mouse RBP4 monoclonal antibody (15  $\mu$ g/ml) at 37°C 1 hour. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-rabbit IgG secondary antibody (green).