

Anti-PCSK9 antibody



Catalog Number: 100582

Product name

Anti-PCSK9 antibody

Immunogen

[Human PCSK9 \(His Tag\) recombinant protein](#)

Specificity

Human PCSK9 / NARC1

Antibody description

Rabbit polyclonal to PCSK9

Preparation

Produced in rabbits immunized with purified, recombinant Human PCSK9 / NARC1 (rh PCSK9 / NARC1; NP_777596.2; Met1-Gln692). PCSK9 / NARC1 specific IgG was purified by Human PCSK9 / NARC1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human PCSK9 / NARC1. The detection limit for Human PCSK9 / NARC1 is < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

Validations



PCSK9 / NARC1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry