

# Anti-PEBP1/RKIP antibody



Catalog Number: 100768

## Product name

Anti-PEBP1/RKIP antibody

## Immunogen

[Human PEBP1/RKIP recombinant protein](#)

## Specificity

Human PEBP1

## Antibody description

Rabbit polyclonal to PEBP1/RKIP

## Preparation

Produced in rabbits immunized with purified, recombinant Human PEBP1 (rh PEBP1; NP\_002558.1; Met1-Lys187). PEBP1 specific IgG was purified by Human PEBP1 affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions


WB: 2-10  $\mu$ g/ml

ELISA: 0.1-0.2  $\mu$ g/mL

This antibody can be used at 0.1-0.2  $\mu$ g/mL with the appropriate secondary reagents to detect Human PEBP1. The detection limit for Human PEBP1 is < 0.039 ng/well.


IP: 0.5-2  $\mu$ g/mg of lysate

## Validations



Lanes	A	B
Items		
Sample (whole cell lysate)	A549	Hela
Sample Volume ( $\mu$ g/lane)	30	30
Gel	13% SDS-PAGE reducing gel	
Recommended Concentration	2-10 $\mu$ g/ml	
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution	
Developed using Odyssey imaging system.		
Explanation	Predicted band size : 21 kDa phenylethanolamine band size : 21 kDa	

PEBP1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



Lanes	A	B
Items		
Sample (whole cell lysate)	HeLa	A549
Sample quantity	0.5 mg	
IP antibody quantity	2 $\mu$ g	
Protein G agarose	15 $\mu$ l of 50% Protein G Agarose	
Gel	13% SDS-PAGE reducing gel	
Primary antibody	PEBP1 antibody at 5 $\mu$ g/ml	
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution.	

PEBP1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation