

Anti-Endoglin/CD105 antibody



Catalog Number: 100871

Product name

Anti-Endoglin/CD105 antibody

WB: 1-2 µg/mL

IHC-P: 0.1-2 µg/mL


Specificity

Human CD105

Validations

Antibody description

Rabbit polyclonal to Endoglin/CD105



Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human CD105, and purified by antigen affinity chromatography.

Lanes	A
Items	
Sample (Recombinant protein)	CD105
Sample Volume (ng/lane)	100
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	1-2 µg/ml
Secondary Antibody	HRP Conjugated Goat anti-rabbit IgG (H+L) at 1:1000 dilution

Formulation

0.2 µm filtered solution in PBS with 5% trehalose


Endoglin / CD105 / ENG Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.



Endoglin / CD105 / ENG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Clonality

Polyclonal

Immunochemical staining of human CD105 in human kidney with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to endothelial cells.

Ig Type

Rabbit IgG


Applications

WB, IHC-P

Dilutions


Anti-Endoglin/CD105 antibody

Catalog Number: 100871



Endoglin / CD105 / ENG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human CD105 in human lung with rabbit polyclonal antibody (1 $\mu\text{g}/\text{mL}$, formalin-fixed paraffin embedded sections). Positive staining was localized to endothelial cells.



Endoglin / CD105 / ENG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human CD105 in human tonsil with rabbit polyclonal antibody (1 $\mu\text{g}/\text{mL}$, formalin-fixed paraffin embedded sections). Positive staining was localized to endothelial cells.