

Anti-IFNAR1 antibody



Catalog Number: 100928

Product name

Anti-IFNAR1 antibody

Immunogen

[Human IFNAR1 \(His Tag\) recombinant protein](#)

Specificity

Human IFNAR1

Antibody description

Mouse monoclonal to IFNAR1

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IFNAR1 (rh IFNAR1; P17181-1; Met1-Lys436) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze
!Sodium azide is toxic to cells and should be

disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



IFNAR1 Antibody (PE), Mouse MAb, Flow cytometric analysis

Flow cytometric analysis of Human IFNAR1 expression on Jurkat cells.