

# Anti-PGA4/Pepsinogen A antibody



Catalog Number: 100937

## Product name

Anti-PGA4/Pepsinogen A antibody

## Immunogen

[Human PGA4/Pepsinogen A \(His Tag\) recombinant protein](#)

## Specificity

Human PGA4 / Pepsinogen A

## Antibody description

Rabbit polyclonal to PGA4/Pepsinogen A

## Preparation

Produced in rabbits immunized with purified, recombinant Human PGA4 / Pepsinogen A (rh PGA4 / Pepsinogen A; P00790; Met1-Ala388). PGA4 / Pepsinogen A specific IgG was purified by Human PGA4 / Pepsinogen A affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human PGA4 / Pepsinogen A. The detection limit for Human PGA4 / Pepsinogen A is approximately < 0.039 ng/well.

IHC-P: 0.1-2 µg/ml

## Validations



PGA4 / Pepsinogen A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human PGA in human gastric cancer with rabbit polyclonal antibody (0.2 µg/mL, formalin-fixed paraffin embedded sections).