

Anti-CD132 / IL2RG antibody



Catalog Number: 100941

Product name

Anti-CD132 / IL2RG antibody

Immunogen

[Mouse CD132 / IL2RG \(His Tag\) recombinant protein](#)

Specificity

Mouse CD132 / IL2RG

Antibody description

Rabbit monoclonal to CD132 / IL2RG

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD132 / IL2RG (rM CD132 / IL2RG; NP_038591.1; Met1-Ala263) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD132 / IL2RG Antibody (PE), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of mouse IL2RG(CD132) expression on mouse splenocytes. BABL/c splenocytes were stained with PE-conjugated anti-Mouse IL2RG(CD132). The histogram were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.