

# Anti-Syndecan-4/SDC4 antibody



Catalog Number: 100953

## Product name

Anti-Syndecan-4/SDC4 antibody

## Immunogen

[Mouse Syndecan-4/SDC4 \(His Tag\) recombinant protein](#)

## Specificity

Mouse Syndecan-4 / SDC4

## Antibody description

Rabbit monoclonal to Syndecan-4/SDC4

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Syndecan-4 / SDC4 (rM Syndecan-4 / SDC4; O35988; Met1-Val146) and conjugated with PE under optimum conditions, the unreacted PE was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



Syndecan-4 / SDC4 Antibody (PE), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of mouse SDC4 expression on mouse splenocytes. BABL/c splenocytes were stained with PE-conjugated anti-Mouse SDC4 and FITC-conjugated anti-Mouse CD3 (BD Pharmingen™ Cat. No.553061, left) or APC-conjugated anti-Mouse CD19 (BD Pharmingen™ Cat. No.550992, right).