

# Anti-ADAM10 antibody



Catalog Number: 101003

## Product name

Anti-ADAM10 antibody

## Specificity

Human ADAM10 / CD156c

## Antibody description

Rabbit polyclonal to ADAM10

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human ADAM10 / CD156c, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 10-50 µg/ml

## Validations



Items	Lanes	A
Sample (whole cell lysate)		jurkat
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-50 µg/ml
Secondary Antibody		DyLight 900-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system		
Explanation		Predicted band size: 84 kDa Observed band size: 84 kDa Additional bands at: 50 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).

ADAM10 / CD156c Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot