

Anti-RNase A / Ribonuclease A / RNASE1 antibody



Catalog Number: 101061

Product name

Anti-RNase A / Ribonuclease A / RNASE1 antibody

Immunogen

[Human RNase A / Ribonuclease A / RNASE1 \(His Tag\) recombinant protein](#)

Specificity

Human RNase A / Ribonuclease A / RNASE1

Antibody description

Rabbit polyclonal to RNase A / Ribonuclease A / RNASE1

Preparation

Produced in rabbits immunized with purified, recombinant Human RNase A / Ribonuclease A / RNASE1 (rh RNase A / Ribonuclease A / RNASE1; P07998; Met1-Thr156). RNase A / Ribonuclease A / RNASE1 specific IgG was purified by Human RNase A / Ribonuclease A / RNASE1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 µg/ml

Validations



RNase A / Ribonuclease A / RNASE1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry