

Anti-PAX3 antibody



Catalog Number: 101085

Product name

Anti-PAX3 antibody

Specificity

Human PAX3

Antibody description

Rabbit polyclonal to PAX3

Preparation

Produced in rabbits immunized with purified, a synthetic peptide corresponding to the C-terminus of the human PAX3. Human PAX3 specific IgG was purified by Human PAX3 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Items	Lanes	A	B	C
Sample (whole cell lysate)		Jurkat	A431	A549
Sample Volume (µg/lane)		30	30	30
Gel	13% SDS-PAGE reducing gel			
Recommended concentration		5-10 µg/ml		
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Explanation		Developed using Odyssey imaging system. Predicted band size : 53 kDa Observed band size : 60 kDa		

PAX3 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot