

# Anti-RAB5A antibody



Catalog Number: 101182

## Product name

Anti-RAB5A antibody

## Specificity

human RAB5A

## Antibody description

Rabbit polyclonal to RAB5A

## Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the C-terminus of the human RAB5A, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-2 µg/mL

## Validations



RAB5A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human RAB5A in human small intestine (from 2 donors) with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).