

Anti-RAD18 antibody



Catalog Number: 101212

Product name

Anti-RAD18 antibody

Specificity

Human RAD18

Antibody description

Rabbit polyclonal to RAD18

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human RAD18, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC-P, IF, ICC/IF

Dilutions

IHC-P: 0.1-2 µg/mL


ICC/IF: 0.5-1.5 µg/mL

Validations



RAD18 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human RAD18 in human testis with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded



sections).

RAD18 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence

Immunofluorescence staining of RAD18 in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-human RAD18 polyclonal antibody (1 µg/ml) at 4°C overnight. Then cells were