

# Anti-ADAM8(CD156) antibody



Catalog Number: 101369

## Product name

Anti-ADAM8(CD156) antibody

## Specificity

Human CD156a

## Antibody description

Rabbit polyclonal to ADAM8(CD156)

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the human CD156a, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 0.1-2 µg/ml

## Validations



CD156a Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human ADAM8 in human placenta with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded



sections).

CD156a Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human ADAM8 in human kidney with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).