

Anti-Coagulation Factor III / Tissue Factor / CD142 antibody



Catalog Number: 101386

Product name

Anti-Coagulation Factor III / Tissue Factor / CD142 antibody

Immunogen

[Mouse Coagulation Factor III / Tissue Factor / CD142 \(His Tag\) recombinant protein](#)

Specificity

Mouse Coagulation Factor III / Tissue Factor / CD142

Antibody description

Rabbit monoclonal to Coagulation Factor III / Tissue Factor / CD142

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Coagulation Factor III / Tissue Factor / CD142 (rM Coagulation Factor III / Tissue Factor / CD142; NP_034301.3; Met1-Glu251).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 1.5-3 µg/Test

Validations



Coagulation Factor III / Tissue Factor / CD142 Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse F3(CD142) expression on Raw264.7 cells. Cells were stained with purified anti-Mouse F3(CD142), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.