

Anti-HA antibody



Catalog Number: 101457

Product name

Anti-HA antibody

Immunogen

[H1N1 HA \(His Tag\) recombinant protein](#)

Specificity

Influenza A H1N1 (A/California/04/2009)
Hemagglutinin / HA

Antibody description

Rabbit polyclonal to HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA (rh Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA; ACP41105.1; Met1-Gln529). Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA specific IgG was purified by Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of

receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, IP

Dilutions

WB: 0.5-1 µg/ml

Validations



Items	Lanes	A	B
Sample (Protein)		H1N1-C-Hah	H1N1-C-Hah
Sample Volume (µg/lane)		100	10
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-1 µg/ml		
Secondary Antibody	Goat anti Rabbit IgG (H+L)/HRP:0.4µg/ml.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 75-80 kDa Observed band size : 75-80 kDa		

Influenza A H1N1 (A/California/04/2009) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot