

# Anti-PAX7/HUP1 antibody



Catalog Number: 101615

## Product name

Anti-PAX7/HUP1 antibody

## Specificity

Human PAX7 / HUP1

## Antibody description

Rabbit polyclonal to PAX7/HUP1

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human PAX7 / HUP1, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-20 µg/ml

## Validations



Items	Lanes	A
Sample (whole cell lysate)		Hela
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-20 µg/ml
Secondary Antibody		Dylight 800 Labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
		Developed using Odyssey imaging system.
Explanation		Predicted band size : 57 kDa Observed band size : 55 kDa Additional bands at : 35 kDa and 32 kDa (We are unsure as to the identity of these extra bands.)

PAX7 / HUP1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot