

Anti-RELT/TNFRSF19L antibody



Catalog Number: 102140

Product name

Anti-RELT/TNFRSF19L antibody

Immunogen

[Human RELT/TNFRSF19L \(His Tag\) recombinant protein](#)

Specificity

Human RELT / TNFRSF19L

Antibody description

Mouse monoclonal to RELT/TNFRSF19L

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human RELT / TNFRSF19L (rh RELT / TNFRSF19L; NP_116260.2; Met1-Ala160). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

IF, ICC/IF

Dilutions

ICC/IF: 10-25 µg/mL

Validations



RELT / TNFRSF19L Antibody, Mouse MAb, Immunofluorescence

Confocal immunofluorescence analysis of Human RELT in SKBR3 cells. Cells were fixed with 4% PFA, permeabilized with 1% Triton X-100 in PBS, blocked with 10% serum, and incubated with Mouse anti-Human RELT monoclonal antibody (15 µg/ml) at 37°C 1 hour. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-mouse IgG secondary antibody (green). Positive staining was localized to cytoplasm and plasma membrane.