

Anti-PILR-alpha/PILRA antibody



Catalog Number: 102888

Product name

Anti-PILR-alpha/PILRA antibody

Immunogen

[Mouse PILR-alpha/PILRA \(His Tag\) recombinant protein](#)

Specificity

Mouse PILRA

Antibody description

Rabbit polyclonal to PILR-alpha/PILRA

Preparation

Produced in rabbits immunized with purified, recombinant Mouse PILRA (rM PILRA; NP_705730.1; Met 1-Val 197). PILRA specific IgG was purified by mouse PILRA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse PILRA. The detection limit for Mouse PILRA is approximately 0.00245 ng/well.

IHC-P: 0.1-1 µg/mL

Validations



PILR-alpha / PILRA Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse PILRA in mouse intestine with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).