

Anti-CD20/MS4A1 antibody



Catalog Number: 103691

Product name

Anti-CD20/MS4A1 antibody

Immunogen

[Ferret CD20/MS4A1 \(His Tag\) recombinant protein](#)

Specificity

Ferret CD20 / MS4A1

Antibody description

Rabbit polyclonal to CD20/MS4A1

Preparation

Produced in rabbits immunized with purified, recombinant Ferret CD20 / MS4A1 (rF CD20 / MS4A1; Glu 213-Pro 297). CD20 / MS4A1 specific IgG was purified by Ferret CD20 / MS4A1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Ferret CD20. The detection limit for Ferret CD20 is approximately 0.00245 ng/well.


IHC-P: 0.1-1 µg/mL

Validations



CD20 / MS4A1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of ferret CD20 in ferret spleen with rabbit polyclonal antibody (0.1 µg/mL, formalin-fixed paraffin embedded sections).



CD20 / MS4A1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Anti-CD20/MS4A1 antibody



Catalog Number: 103691

Immunochemical staining of ferret CD20 in ferret colon with rabbit polyclonal antibody (0.1 µg/mL, formalin-fixed paraffin embedded sections).