

Anti-CD200 antibody



Catalog Number: 103828

Product name

Anti-CD200 antibody

Specificity

Mouse CD200

Antibody description

Rabbit monoclonal to CD200

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD200 (rM CD200; NP_001076017.1; Met 1-Gly 232).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 0.5-2 µg/Test

Validations



CD200 Antibody, Rabbit MAb, Flow cytometric

Flow cytometric analysis of CD200 expression on spleen lymphocytes. BALB/c splenocytes were stained with Mouse CD200 Monoclonal Ab (50074-R002, 1 µg/test), then a FITC-conjugated second step antibody. The histogram were derived from the gated events based on light scattering characteristics of lymphocytes.