

# Anti-STX7 / Syntaxin 7 antibody



Catalog Number: 104202

## Product name

Anti-STX7 / Syntaxin 7 antibody

## Immunogen

[Rat STX7 / Syntaxin 7 \(His Tag\) recombinant protein](#)

## Specificity

Rat STX7 / Syntaxin 7

## Antibody description

Rabbit polyclonal to STX7 / Syntaxin 7

## Preparation

Produced in rabbits immunized with purified, recombinant Rat STX7 / Syntaxin 7 (rR STX7 / Syntaxin 7; NP\_068641.2; Ser2-Leu238). STX7 / Syntaxin 7 specific IgG was purified by Rat STX7 / Syntaxin 7 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Rat STX7 / Syntaxin 7.

IHC-P: 0.1-2 µg/mL

## Validations



STX7 / Syntaxin 7 Antibody, Rabbit PAb, Antigen Affinity Purified

Immunochemical staining of rat STX7 in rat kidney with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).