

Anti-IL-6 antibody



Catalog Number: 104410

Product name

Anti-IL-6 antibody

Immunogen

[Mouse IL-6 recombinant protein](#)

Specificity

Mouse IL6 / IL-6

Antibody description

Rabbit polyclonal to IL-6

Preparation

Produced in rabbits immunized with purified, recombinant Mouse IL6 / IL-6 (rh IL6 / IL-6; P08505; Phe25-Thr211). IL6 / IL-6 specific IgG was purified by Mouse IL6 / IL-6 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-50 µg/mL

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse IL6 / IL-6. The detection limit for Mouse IL6 / IL-6 is < 0.039 ng/well.

Validations



IL6 / IL-6 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot