

Anti-Semaphorin 4D/SEMA4D/CD100 antibody



Catalog Number: 104497

Product name

Anti-Semaphorin 4D/SEMA4D/CD100 antibody

Immunogen

[Human Semaphorin 4D/SEMA4D/CD100 \(His Tag\) recombinant protein](#)

Specificity

Human SEMA4D / CD100

Antibody description

Mouse monoclonal to Semaphorin 4D/SEMA4D/CD100

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SEMA4D / CD100 (rh SEMA4D / CD100; Q92854-1; Met1-Arg734) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Semaphorin 4D / SEMA4D / CD100 Antibody (PE), Mouse MAb, Flow cytometric

Profile of peripheral blood lymphocytes analyzed by flow cytometry.