

Anti-IL16 antibody

Catalog Number: 104579

Product name

Anti-IL16 antibody

Immunogen

[Human IL16 \(His Tag\) recombinant protein](#)

Specificity

Human IL-16

Antibody description

Rabbit monoclonal to IL16

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human IL-16 (rh IL-16; AAC12732.1; Pro2-Ser130).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



IL-16 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human IL16 in human spleen with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

IL-16 Antibody, Rabbit MAb,
Immunohistochemistry

Immunochemical staining of human IL16 in human tonsil with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).