

Anti-CD200 antibody



Catalog Number: 104586

Product name

Anti-CD200 antibody

Immunogen

[Mouse CD200 \(His Tag\) recombinant protein](#)

Specificity

Mouse CD200

Antibody description

Rabbit monoclonal to CD200

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD200 (rM CD200; NP_001076017.1; Met1-Gly232) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be

disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD200 Antibody (APC), Rabbit MAb, Flow cytometric analysis

Analysis of CD200 expression on spleen lymphocytes.