

Anti-DDX4 / VASA antibody



Catalog Number: 104612

Product name

Anti-DDX4 / VASA antibody

Specificity

Human DDX4

Antibody description

Mouse monoclonal to DDX4 / VASA

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with a synthetic peptide corresponding to the C-terminus of the human DDX4. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

IHC-P

Dilutions

IHC-P: 10-30 µg/mL

Validations



DDX4 Antibody, Mouse MAb,
Immunohistochemistry

Immunochemical staining of human DDX4 in human testis with mouse monoclonal antibody (15 µg/mL, formalin-fixed paraffin embedded sections).