

# Anti-HA antibody



Catalog Number: 104858

## Product name

Anti-HA antibody

## Immunogen

[H1N1 HA \(His Tag\) recombinant protein](#)

## Specificity

Influenza A H1N1 (A/Puerto Rico/8/34)  
Hemagglutinin / HA

## Antibody description

Rabbit polyclonal to HA

## Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/Puerto Rico/8/34) Hemagglutinin / HA (ABD77675.1; Met1-Gln528). Influenza A H1N1 (A/Puerto Rico/8/34) Hemagglutinin / HA specific IgG was purified by Influenza A H1N1 (A/Puerto Rico/8/34) Hemagglutinin / HA affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC-P


## Dilutions

WB: 0.5-1 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza A H1N1 (A/Puerto Rico/8/34) Hemagglutinin / HA.

## Validations



Lanes	A	B	C
Items			
Sample (Recombinant protein)	H1N1-PR-34-HA-His (Cat#11684-V08H)		
Sample Volume (ng/lane)	20	5	1
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-1 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.		
Developed using Odyssey imaging system.			

Influenza A H1N1 (A/Puerto Rico/8/34) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot