

Anti-Androgen receptor/NR3C4 antibody



Catalog Number: 104861

Product name

Anti-Androgen receptor/NR3C4 antibody

Specificity

Human Androgenreceptor / AR

Antibody description

Rabbit polyclonal to Androgen receptor/NR3C4

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human Androgenreceptor / AR, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Lanes	A	B	C
Items			
Sample (whole cell lysate)	HepG2	Hela	293T
Sample Volume (µg/lane)	30	30	30
Cell	7.5% SDS-PAGE reducing gel		
Recommended Concentration	5-10 µg/ml		
Secondary Antibody	Dylight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
	Developed using Odyssey imaging system.		
Explanation	Predicted band size : 98 kDa Observed band size : 100 kDa Additional bands at : 82 kDa [We are unsure as to the identity of these extra bands.]		

Androgenreceptor / AR Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot