

Anti-NFKBIB / IκB-beta antibody



Catalog Number: 104912

Product name

Anti-NFKBIB / IκB-beta antibody

Specificity

Human NFKBIB / IκB-beta

Antibody description

Rabbit polyclonal to NFKBIB / IκB-beta

Preparation

Produced in rabbits immunized with purified, recombinant Human NFKBIB / IκB-beta (rh NFKBIB / IκB-beta; AAH15528.1; 1M-356V). NFKBIB / IκB-beta specific IgG was purified by Human NFKBIB / IκB-beta affinity chromatography.

Formulation

0.2 μm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP


Dilutions

ELISA: 0.1-0.2 μg/ml

This antibody can be used at 0.1-0.2 μg/ml with the appropriate secondary reagents to detect Human NFKBIB / IκB-beta.

IP: 0.5-2 uL/mg of lysate

Validations



Items	Label	A
Sample (whole cell lysate)		Jurkat
Sample quantity		0.5 mg
IP antibody quantity		2 μg
Immunomagnetic beads Protein G		60 μg
Gel		13% SDS-PAGE reducing gel
Primary antibody		GST-NFKBIB antibody at 10 μg/ml
Secondary antibody		Cleav-Host IP Detection Reagent (HRP) at 1:500 dilution

NFKBIB / IκB-beta Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation