

# Anti-KAZALD1 antibody



Catalog Number: 104931

## Product name

Anti-KAZALD1 antibody

## Specificity

Human KAZALD1

## Antibody description

Rabbit polyclonal to KAZALD1

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human KAZALD1, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-10 µg/ml

## Validations



Lanes	A	B	C	D
Items				
Sample (whole cell lysate)	HepG2	Hela	A431	THP1
Sample Volume (µg/lane)	30	30	30	30
Get	13% SDS-PAGE reducing gel			
Recommended Concentration	5-10 µg/ml			
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.			
	Developed using Odyssey imaging system.			
Explanation	Predicted band size : 33kDa Observed band size : 45 kDa Additional bands at : 39 kDa and 37 kDa (We are unsure as to the identity of these extra bands.)			

KAZALD1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot