

Anti-ADAM17 antibody



Catalog Number: 104968

Product name

Anti-ADAM17 antibody

Specificity

Human ADAM17 / CD156b

Antibody description

Rabbit polyclonal to ADAM17

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human CD156b, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 4-6 µg/mg of lysate

Validations



Item	Lanes	A
Sample (whole cell lysate)		Raw2647
Sample quantity		0.5 mg
IP antibody quantity		2 µg
Protein G agarose		15 µl of 50% Protein G Agarose
Cell		13% SDS-PAGE reducing gel
Primary antibody		KLJF-sinoA403 antibody at 5 µg/ml
Secondary antibody		DyLight 800-labeled antibody to rabbit IgG (2H12), at 1:5000 dilution.

ADAM17 / CD156b Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation