

# Anti-PMS2 antibody



Catalog Number: 104975

## Product name

Anti-PMS2 antibody

## Specificity

Human PMS2

## Antibody description

Rabbit polyclonal to PMS2

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human PMS2, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-10 µg/ml

## Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HepG2	HeLa
Sample Volume: (µg/lane)		30	30
Gel		7.5% SDS-PAGE reducing gel	
Recommended Concentration		5-10 µg/ml	
Secondary Antibody		Dylight 800 labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.	
Explanation		Developed using Odyssey imaging system. Predicted band size : 96 kDa Observed band size : 120 kDa Additional bands at : 67 kDa [We are unsure as to the identity of these extra bands. They are possible non-specific bindings].	

PMS2 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot