

Anti-IL-6R/CD126 antibody



Catalog Number: 105053

Product name

Anti-IL-6R/CD126 antibody

Immunogen

[Rat IL-6R/CD126 \(His Tag\) recombinant protein](#)

Specificity

Rat IL-6R / CD126

Antibody description

Rabbit polyclonal to IL-6R/CD126

Preparation

Produced in rabbits immunized with purified, recombinant Rat IL-6R / CD126 (rR IL-6R / CD126; NP_058716.2; Met1-Pro364). IL-6R / CD126 specific IgG was purified by Rat IL-6R / CD126 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-20 µg/ml

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat IL-6R / CD126. The detection limit for Rat IL-6R / CD126 is < 0.039 ng/well.

Validations



IL-6R / CD126 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot