

# Anti-caspase-3 / CASP-3 antibody



Catalog Number: 105079

## Product name

Anti-caspase-3 / CASP-3 antibody

## Specificity

Human CASP3

## Antibody description

Rabbit polyclonal to caspase-3 / CASP-3

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human CASP3, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB

## Dilutions

WB: 5-10 µg/ml

IP: 4-6 µg/mg of lysate

## Validations



Items	Lanes	A
Sample (whole cell lysate)		NIH3T3
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-10 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size : 32 kDa Observed band size : 32 kDa

CASP3 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot