

Anti-G-CSFR/CD114 antibody



Catalog Number: 105302

Product name

Anti-G-CSFR/CD114 antibody

Specificity

Human G-CSFR / CD114 / CSF3R

Antibody description

Rabbit monoclonal to G-CSFR/CD114

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human G-CSFR / CD114 / CSF3R (rh VE-Cadherin / CD144 / CDH5; CAA56306.1; Met1-Gln593) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



G-CSFR / CD114 / CSF3R Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Human CD144 expression on HUVEC cells. Cells were stained with FITC-conjugated anti-Human CD144. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.