

# Anti-Vimentin antibody



Catalog Number: 105334

## Product name

Anti-Vimentin antibody

## Specificity

Human VIM

## Antibody description

Rabbit polyclonal to Vimentin

## Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the C-terminus of the human Vimentin., and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions


IHC-P: 0.1-2 µg/mL

## Validations



Vimentin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human Vimentin in human malignant melanoma with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).



Vimentin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human Vimentin in human colon carcinoma with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).