

Anti-IL-9 antibody

Catalog Number: 105565

Product name

Anti-IL-9 antibody

Specificity

Mouse IL9

Antibody description

Rabbit polyclonal to IL-9

Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the mouse IL9 , and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 10-20 µg/mL


IHC-P: 0.1-2 µg/mL

Validations



IL9 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse IL9 in mouse spleen (from 2 donors) with rabbit polyclonal antibody (0.2 µg/mL, formalin-fixed paraffin



embedded sections).

IL9 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot