

Anti-Prostasin/PRSS8 antibody



Catalog Number: 105985

Product name

Anti-Prostasin/PRSS8 antibody

Immunogen

[Human Prostasin/PRSS8 \(His Tag\) recombinant protein](#)

Specificity

Human Prostasin / PRSS8

Antibody description

Rabbit monoclonal to Prostasin/PRSS8

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Prostasin / PRSS8 (rh Prostasin / PRSS8; NP_002764.1; Ala30-Arg322).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 2-10 µg/mL

Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (µg/lane)		30
Cell		13% SDS-PAGE reducing gel
Recommended Concentration		2.5 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size: 38 kDa Observed band size: 38 kDa

Prostasin / PRSS8 Antibody, Rabbit MAb, Western blot