

# Anti-HPGD/15-PGDH antibody



Catalog Number: 106039

## Product name

Anti-HPGD/15-PGDH antibody

## Immunogen

[Mouse HPGD/15-PGDH \(His Tag\) recombinant protein](#)

## Specificity

Mouse HPGD / 15-PGDH

## Antibody description

Rabbit monoclonal to HPGD/15-PGDH

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse HPGD / 15-PGDH (rM HPGD / 15-PGDH; Q8VCC1; Met1-Ser269).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions


WB: 2-10 µg/mL

IP: 1-4 µg/mg of lysate

## Validations



## HPGD / 15-PGDH Antibody, Rabbit MAb, Western



blot

## HPGD / 15-PGDH Antibody, Rabbit MAb, Immunoprecipitation